Nonwoven Medical Tape	
Introduction	Nonwoven Medical Tape Non e is based on Tan Rayon Fiber Nonwoven. Coated with non-tackified Acrylic-based, developed for skin contact,
Application	Fixing gauze, wound dressings, and bandages, fixing IV lines, catheters, or monitoring equipment, post-operative care and mild wound management, for light fixation and support for medical tape applications.
Property	
Product Code	P9101
Technical Parameters	
Color	brown
Total Thickness[mm]	0.14±0.02
Backing	Tan Rayon Fiber Nonwoven
Adhesive	Non-tackified Acrylic-based, developed for skin contact
180° Peel Strength[N/25]	6. 6
Moisture Vapor Transmission Rate[hrs]	24
Tensile Strength, MD, at break[N/25.4 mm]	22. 2
Storage and shelf life	It is recommended to store in a clean and dry place, away from direct sunlight to prevent damage to the packaging, and to avoid storing with volatile solvents. Store at 10-27° C and humidity at 40~60%.